



Rehabilitation Manufacturing Services Ltd  
 Thompson House  
 Unit 10, Styles Close, Sittingbourne, Kent ME10 3BF  
 Tel: 01795 477280 Fax: 01795 229692  
 Email: sales@rms-kent.co.uk web: www.ineedawheelchair.co.uk

## Doccia Shower Chair

January 2018

Client Name:

Date:

Contact Name:

Tel number:

Address:

### Standard Specification

- Aluminium frame with a non-toxic white epoxy varnish finish
- Crystal coloured washable PVC-elastic straps adapts to the users shape
- Padded headrest separate from the backrest, angle and height adjustable
- 30 ° backrest angle adjustment
- Seat with a 15° tilt, in order to prevent the user from sliding forward
- Detachable and turning sideways sides, in order to facilitate the user's transfers
- Leg rest adjustable in tilt, one piece flip up footrest
- 12 cm puncture proof wheels, rear wheels braked
- Weight capacity 18Kg
- Maximum user weight 110Kg



Dimensions	Seat depth	Seat width	Backrest height	Backrest width	Seat / footrest height	Total width	Total length	Total height
<b>Sizes min./max</b>	40/44 cm	45 cm	75/100 cm	60 cm	36/46 cm	49 cm	92/108 cm	117/145 cm

Part No.	Description	
107104E	203-5 Shower chair	<input type="checkbox"/>

