

# Otto Bock®

QUALITY FOR LIFE



## Skippi

### Experience the World

Exceed your limits, increase your independence, conquer the outdoors: Skippi makes it fun to discover something new every day. Skippi's speed and high degree of maneuverability help it support children's urge to move and the pleasure they derive from making discoveries. Mobility expands your horizons. Children learn by being actively involved in their surroundings. Arm and back rest adjustments are variable. In this way, the chair can be set to facilitate „seeing and doing“. And this goes for traveling too: Skippi is easy to disassemble and fits in any small car. What's more: Skippi grows with the child! Adapting its seat depth and width is child's play.

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt · Phone +49 5527 848-1590 · Fax +49 5527 848-1676 · [healthcare@ottobock.de](mailto:healthcare@ottobock.de) · [www.ottobock.com](http://www.ottobock.com)



# Skippy – Power Wheelchair for children

## Options

- Seat tilt mechanism, electric or with gas compression spring
- electrical Backrest angle adjustment
- Seating shell adapter
- Tray
- Control for attendant
- Bumber bar
- Single-panel footrest with aluminum footplate

## Technical Data

- Max. speed .....6 km/h
- Battery capacity ..... 12V, 31Ah C20
- Distance range ..... approx. 30 km
- Max. load capacity .....50 kg
- Empty weight .....68 kg
- Max. obstacle height that can be cleared ...50 mm
- Overall width .....57 cm
- Overall length without legrests .....70 cm
- Overall length with legrests .....85 cm
- Turning radius.....55 cm
- Seat height .....43 cm
- Seat width .....26 – 38 cm
- Seat depth.....30 – 37 cm
- Seat inclination (optional: continuous).....-9°/3°/12° (0°– 25°)
- Armrest height .....16 – 27 cm
- Back height .....35/40/42/44/46 cm
- Back width .....30 cm
- Back angle (optional: continuous) .....0-25°
- Lower leg length .....11 – 26 cm
- Frame colors .....yellow, blue, red



Height adjustable push handles for attendant



Electrical backrest angle adjustment



Seat adaptable to size



Detachable footrests

Supplied By:



Medway House, 277 Gillingham Road,  
Gillingham, Kent, ME7 4QX  
Tel: 01634 578881 Fax: 01634 581456  
Email: sales@rms-kent.co.uk