



Thompson House, Unit 10, Styles Close, Sittingbourne, Kent, ME10 3BF.

Tel: +44 (0) 1795 477280 **Fax:** +44 (0) 1795 229692 Email: sales@rms-rehab.co.uk

www.rms-rehab.co.uk



sales@rms-rehab.co.uk 01795 477280

www.rms-rehab.co.uk

STANDSURE

SIT-TO-STAND THERAPY AID





SIT-TO-STAND THERAPY AID



The first ever therapy tool specifically designed to be incorporated into rehabilitation programs for users with difficulties maintaining the neutral balance of their body while walking or standing, minimising the risk of accidents while carrying out therapy.

Using StandSure means the user experiences a greater level of comfort and stability while practicing a repetitive sit-to-stand motion, keeping the patients feet safely grounded throughout the session.

Manufactured in three sizes suitable for both children and adults - Standsure is straightforward, lightweight, easy to carry and easy to store - whilst showing incredible results in posture improvement, balance, range of motion, and upper body limb function.

KEY FEATURES

- > Available in three varying sizes suitable for both adults and children
- > Made from extremely durable materials that are lightweight and easy to carry
- > An affordable and easy-to-use therapy aid
- > Includes ankle cuffs + interchangeable toe straps
- > Easy to disinfect for multiple users
- > Very portable for easy storage almost anywhere!



EASY AND SIMPLE TO USE



ERGONOMIC THERAPY AID

Standsure is manufactured to adhere securely to practically any surface. It fits neatly between the legs of chairs or other medical devices such as walkers, standers and therapy room furniture.

EASY AND QUICK SET UP

The board contains width adjustable straps that can be easily and quickly set up either by a healthcare professional, career, or the user - when directed to be used as a stretching routine at home.

LIGHTWEIGHT AND PORTABLE

The StandSure board is lightweight and portable for healthcare professionals to easily carry and use anywhere - perfect for commuting when you need to take sessions outside of the clinic room!

RMS ANKLE STRAPS

NON-SLIP BASE

Standsure has been tested and approved to be a robust and long-lasting product for users of any age and weight, with a base that can be adhered securely to almost any work surface.

01795 477280

-

www.rms-rehab.co.uk





The aid comes with adjustable ankle straps to help the user keep their feet grounded while practising repetitive sit-to-stand exercises, offering comfort and security for both the user and professional

THERAPY BENEFITS

- Stabilises the user by grounding the feet in a healthy position while practising repetitive sit-to-stand or any related rehabilitation program.
- When included in therapy sessions, StandSure can be used as an 'extra pair of hands' for therapists who work alone.
- A comfortable therapy aid designed for users of all ages.
- Repetitive sit-to-stand helps to improve posture, balance, range of movement, upper body function, gait, flexibility and endurance.
- Reduces asymmetrical weight bearing, benefiting the whole-body system, potentially helping to reduce the risk of falls.
- If used early enough, can help prevent the development of contracting issues.
- Aid professionally and technically approved.
- Product is fully safety tested under ISO 9004:2015.

TECHNICAL SPECIFICATIONS

	SMALL	MEDIUM	LARGE
OVERALL BOARD SIZE	45 cm x 45 cm	52.5 cm x 52.5 cm	60 cm x 60 cm
DISTANCE BETWEEN FEET	10.5 cm	17.5 cm	24.5 cm
USER FOOT WIDTH	8 cm MAX	9.5 cm MAX	10.5 cm MAX
USER FOOT LENGTH	13 cm to 17 cm	15 cm to 22 cm	15 cm to 23 cm
ANKLE CIRCUMFERENCE	16.5 cm to 21 cm	17.5 cm to 24 cm	24 cm to 28 cm