### SPECIFICATIONS

additional specifications at zingstanders.com

Size 1 - 70 lbs. (32kg)   Size 2 - 154 lbs. (70kg)		
Transfer Height  Transfer Height  Transfer Height  Size 2- 40"-60" (102-152cm)  Size 1- 32" (81cm)  Size 2- 31" (79cm)  Size 2- 31" (79cm)  Size 2- 15"-33" (38-84cm)  Size 2- 15"-33" (38-84cm)  Size 2- 6" (15cm)  Size 1- 5" (13cm)  Size 1- 5" (13cm)  Size 1- 5" (13cm)  Size 2- 7" (18cm)  Weight of Basic Unit  Size 1- 55 lbs (25kg)  Size 2- 82 lbs (37kg)  Size 1- 20"x28" (51x71cm)  Size 2- 26.75"x39.75" (68x101cm)  Available on all models  O°-90°  (not available on Vertical & Prone)  Prone  Prone  Prone  Positions  Positions  Positions  Positions  Positions  Positions  Swing-Away Tray  Dual Control  Gas Spring Lockout  Transfer Aid  Transfer Aid  Tool Free Adjustment  Yes	Weight Limit	=
Hip to Foot Plate Range  Upper Body Support Height Range  Knee Pad Depth Range  Knee Pad Depth Range  Knee Pad Depth Range  Weight of Basic Unit  Frame Footprint  Size 1- 5° (15cm)  Size 1- 5° (15cm)  Size 2- 7° (18cm)  Size 2- 7° (18cm)  Size 2- 8° (51x71cm)  Size 1- 20° x28° (51x71cm)  Size 2- 8° (51x71cm)  Size 2- 8° (15cm)  Size 1- 5° Ibs (25kg)  Size 1- 20° x28° (51x71cm)  Size 2- 8° (15cm)  Size 1- 5° (15cm)  Size 1- 5° Ibs (25kg)  Size 1- 20° x28° (51x71cm)  Size 2- 8° (15cm)  Size 1- 20° (15cm)  Size 1- 20° x28° (15cm)  Size 1- 20° (20° (10cm)  Size 1- 20° (20° (	Height Range	. ,
Upper Body Support Height Range  Knee Pad Depth Range  Knee Pad Depth Range  Knee Pad Depth Range  Weight of Basic Unit  Size 1 - 5" (13cm)  Size 2 - 7" (18cm)  Size 2 - 7" (18cm)  Size 2 - 82 lbs (25kg)  Size 2 - 82 lbs (37kg)  Frame Footprint  30° Independent Leg Abduction on each leg  Supine  Prone  Prone  Prone  O°-90° (not available on Vertical & Prone)  Positions  Positions  Positions  Positions  Adjustable Tray  Adjustable Tray  Available on all models  Supine (size 1 & 2)  Vertical (size 1)  Prone (size 1 & 2)  Vertical (size 1)  Prone (size 1 & 2)  Available on all models  Available on all models  Available on all models  Supine (Size 1 & 2)  Vertical (size 1)  Prone (size 1 & 2)  Vertical (size 1)  Prone (size 1 & 2)  Available on all models  Available on all models  Available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Available (not available on Vertical)  Flat-to-Load (not available on Vertical)  Transfer Aid  Flat-to-Load (not available on Vertical)  Yes	Transfer Height	
Size 2-6" (15cm)	Hip to Foot Plate Range	
Weight of Basic Unit  Weight of Basic Unit  Size 1- 55 lbs (25kg) Size 2- 82 lbs (37kg)  Frame Footprint  Size 1- 20"x28" (51x71cm) Size 2- 26.75"x39.75" (68x101cm)  Available on all models  Supine  O°-90° (not available on Vertical & Prone)  Prone  Prone  Prone  O°-20°-MPS, MPS TT; O°-90°-Prone (not available on Supine, Supine TT, & Vertical)  E0641-Multi-Position Stander MPS, MPS TT  E0638-Single Position Stander MPS, MPS TT  E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Available on all models  Available on all models  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Dual Control  Gas Spring Lockout  Available (not available on Vertical)  Flat-to-Load (not available on Vertical)  Transfer Aid  Flat-to-Load (not available on Vertical)  Yes	Upper Body Support Height Range	, ,
Size 2- 82 lbs (37kg)  Frame Footprint  Size 1- 20"x28" (51x71cm) Size 2- 26.75"x39.75" (68x101cm)  Available on all models  Supine  O°-90° (not available on Vertical & Prone)  Prone  Prone  Prone  HCPCS  Fositions  Positions  Positions  Positions  Prone (size 1 & 2) Vertical (size 1) Prone (size 1 & 2) Vertical (size 1) Prone (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Available on all models  Swing-Away Tray  Available on all models  Swing-Away Tray  Available on all models  Available on all models  Fadded Strapping  Available (not available on Prone and Vertical)  Available (not available on Vertical)  Fadded Strapping  Available (not available on Vertical)  Fadded Strapping  Available (not available on Vertical)  Fadded Strapping  Available (not available on Vertical)  Fat-to-Load (not available on Vertical)  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Knee Pad Depth Range	
Size 2- 26.75"x39.75" (68x101cm)  30° Independent Leg Abduction on each leg  Supine  0°-90° (not available on Vertical & Prone)  Prone  0°-20°-MPS, MPS TT; 0°-90°-Prone (not available on Supine, Supine TT, & Vertical)  E0641-Multi-Position Stander MPS, MPS TT  E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Supine (size 1 & 2)  Vertical (size 1)  Prone (size 1 & 2)  Adjustable Tray  Available on all models  Swing-Away Tray  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Weight of Basic Unit	
Supine  Supine  O°-90° (not available on Vertical & Prone)  Prone  O°-20°-MPS, MPS TT; O°-90°-Prone (not available on Supine, Supine TT, & Vertical)  E0641-Multi-Position Stander MPS, MPS TT  E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Adjustable Tray  Available on all models  Swing-Away Tray  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Dual Control  Available (not available on Vertical)  Fadded Strapping  Available (not available on Vertical)  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Frame Footprint	· , ,
Prone  O°-20°-MPS, MPS TT; O°-90°-Prone (not available on Supine, Supine TT, & Vertical)  E0641-Multi-Position Stander MPS, MPS TT E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Positions  Positions  Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Vertical (size 1) Prone (size 1 & 2)  Adjustable Tray  Available on all models  Swing-Away Tray  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	30° Independent Leg Abduction on each leg	Available on all models
HCPCS  HCPCS  HCPCS  E0641-Multi-Position Stander MPS, MPS TT E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Supine, Supine TT, Vertical, & Prone  Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Available on all models  Available on all models  Swing-Away Tray  Available on Prone and Vertical) Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Supine	
HCPCS  MPS, MPS TT  E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Positions  Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Adjustable Tray  Available on all models  Swing-Away Tray  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Prone	0°-20°-MPS, MPS TT; 0°-90°-Prone (not available on Supine, Supine TT, & Vertical)
E0638-Single Position Stander Supine, Supine TT, Vertical, & Prone  Positions  Supine (size 1 & 2) Vertical (size 1) Prone (size 1 & 2)  Adjustable Tray  Available on all models  Swing-Away Tray  Available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Available (not available on Vertical)  Gas Spring Lockout  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	HCPCS	
Positions  Vertical (size 1) Prone (size 1 & 2)  Adjustable Tray  Available on all models  Swing-Away Tray  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Transfer Aid  Tool Free Adjustment  Vertical (size 1) Prone (size 1 & 2)  Available (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Flat-to-Load (not available on Vertical)  Yes		
Swing-Away Tray  Available (not available on Prone and Vertical)  Support  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Available (not available on Vertical)  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Positions	Vertical (size 1)
Swing-Away Tray  (not available on Prone and Vertical)  Individual Positioning Components, Padded Strapping  Available (not available on Vertical)  Gas Spring Lockout  Available (not available on Vertical)  Transfer Aid  Flat-to-Load (not available on Vertical)  Tool Free Adjustment  Yes	Adjustable Tray	Available on all models
Support Padded Strapping  Dual Control Available (not available on Vertical)  Gas Spring Lockout Available (not available on Vertical)  Transfer Aid Flat-to-Load (not available on Vertical)  Tool Free Adjustment Yes	Swing-Away Tray	
Gas Spring Lockout  Transfer Aid  Tool Free Adjustment  Tool Spring Lockout  Available on Vertical)  Flat-to-Load (not available on Vertical)  Yes	Support	
Transfer Aid (not available on Vertical)  Transfer Aid Flat-to-Load (not available on Vertical)  Tool Free Adjustment Yes	Dual Control	
Tool Free Adjustment  Yes  Yes	Gas Spring Lockout	
·	Transfer Aid	
Pow'r Up Lift Option Size 2 Only	Tool Free Adjustment	Yes
	Pow'r Up Lift Option	Size 2 Only

#### **Full Line of Standers**

• MPS TT (tilt table)

Supine

Prone

· Supine TT

Vertical

#### **Hygienic Covers**

- · Available in 5 fun patterns for the pelvic, upper body, lateral, head and tray areas
- · Prevents stander from getting soiled, especially during feeding

Zing MPS TT-Size 1

· Easily removed and washable





SEA CREATURES







#### **Frame Colors**












Altimate Medical, Inc. 262 West First Street Morton, MN 56270 USA Tel: 800.342.8968 507.697.6393 Fax: 877.342.8968 507.697.6900 info@zingstanders.com zingstanders.com

Visit our website at zingstanders.com for complete specifications on all Zing products and options.





# The Zing Advantage







#### **True Multi-Position Stander**

- The MPS and MPS TT allow supine, vertical and prone positioning without flipping foot plates or pads, moving the tray or transferring out and back in the stander
- · Smooth transition from a flat-to-load transfer position, supine, vertical and 20 degrees prone - all in one motion, thanks to the unique pivoting action of the Zing MPS and MPS TT
- Provides therapeutic opportunities not possible with other multi-position standers - use prone for head control and when fatigue sets in, simply move to a slight supine position to rest but remain in the standing position
- Frame will only go supine and upright with Prone Position Lockout engaged

#### **True Hip Pivot Leg Abduction**

- Pivot points are in line with the hips, promoting proper alignment in the hips, knees and feet
- · Feet maintain the angle with the rest of the legs no wedges needed and no additional adjustment to the knees or foot plates
- Up to 30 degrees of abduction per leg accomplished through a single knob
- Promotes musculoskeletal development for children during critical growth years
- Abduction can be adjusted and used throughout the entire positioning range

#### Plenty of Growth

- Size 1: Fits children from infancy to 44", and up to 70lbs
- Size 2: Fits children from 40"-60", and up to 154lbs
- · Can easily grow with the child and adapt to their changing needs
- Many supports, pads and straps are available in multiple sizes
- Comfortable transfer height of 32" makes transfers more manageable as children grow

#### Remarkably Modular

- Can be configured as a simple stander or fit a child with the most complex needs
- Designed to grow with the child options can be added as needed
- · Great for a multi-user facility or school

## SUPINE TO PRONE









Supine Lying - Flat-to-Load

#### **Multi-Use Trav**

- The MPS, Supine, MPS TT and Supine TT feature a multi-use tray that is angle adjustable and swings away either direction
- Tray can remain attached during transfers
- More activities are possible by swinging tray away during standing therapy

#### **Unique Lift Mechanisms**

- · Dual controls are available allowing angle adjustment by foot pedal or hand control
- · Gas Spring Lockout prevents an inadvertent change of position, especially important when around other children
- The Pow'r Up Lift option on Size 2 brings the user to the standing position with the touch of a button

#### **Accommodating Contractures**

- No special modifications are needed to accommodate children with contractures
- Using the multi-adjustable foot plates, multi-adjustable knee pads and the independent calf pads, it is possible to adjust for contractures in one or both sides

#### **Easy Adjustment**

- Tool-free adjustment of all components
- Every adjustment can be made with a simple turn of a knob or ratchet handle





